

# COUNT VERIFICATION MACHINE

PCVM



Engineering Solutions For Healthcare



[ParleGlobalUSA.com](http://ParleGlobalUSA.com)



[DJA-Pharma.com](http://DJA-Pharma.com)

**Daiichi Jitsugyo (America), Inc.**

Address : 939 AEC Drive, Wood Dale, IL 60191

Phone : 630-776-2459

Email : [Rnelson@djausa.com](mailto:Rnelson@djausa.com)

Web : [www.dja-pharma.com](http://www.dja-pharma.com)

## COUNT VERIFICATION MACHINE

PARAMETER	SPECIFICATION
Machine	COUNT VERIFICATION MACHINE
Model	PCVM
Material	SS 316L / Delrin for contact part
Required Power (Approx.)	700 Watt
No. of Track	1 Track
Hopper Capacity	1 liter (2-5 ltr. basin can be accommodated as per requirement)
Electric Spec.	230V, 50Hz, 1PH, 3A
Machine Dimension	300 mm x 510 mm x 780 mm
Weight (Approx.)	50 kg

## Salient Features

- » Touch screen HMI with complete PLC control.
  - » Easy to operate with quick change over from product to product.
  - » Compact module for table top count verification for R&D application and IPQC requirement.
  - » Single track vibrator feeder resembling similar process, to that of the production module type counter.
  - » Anodized aluminum channel for smooth unidirectional transfer of the tablet from the hopper to the counting zone.
  - » Highly accurate counting system for counting & recording the pre-set target.
  - » It has got a Dot Matrix Printer attached with refill roll for batch code, product information & reports.
- » Replaceable funnels as change part, to suite different mouth opening of container.
  - » No bottle no fill arrangement.



Single Track Counter with Hopper



Sensor for Container



HMI Panel



Dot Matrix Printer



[ParleGlobalUSA.com](http://ParleGlobalUSA.com)



[DJA-Pharma.com](http://DJA-Pharma.com)

## Daiichi Jitsugyo (America), Inc.

Address : 939 AEC Drive, Wood Dale, IL 60191

Phone : 630-776-2459

Email : [Rnelson@djausa.com](mailto:Rnelson@djausa.com)

Web : [www.dja-pharma.com](http://www.dja-pharma.com)